

Design of Exposure Systems and Bioagent Inactivation Devices: Aerosol Generation, Measurement and Characterization

Dr. Pratim Biswas
Washington University in St. Louis
pratim.biswas@wustl.edu

**Collaborators: WUSTL: D. Chen, M.H. Lee
SLU: M. Buller, R. Lewis, and J. Schriewer**


 Washington University in St. Louis

AA QRL - Aerosol and Air Quality Research Laboratory

Outline of Presentation


- **Aerosol Generation**
- **Aerosol Measurement for bacteria and viruses; protein molecules and fragments.**
- **Characterization of Exposure Chambers - are the aerosols well distributed, flow modeling, transport in respiratory system and deposition**
- **Inactivation of Bioagents in Airborne Streams - just brief mention.**

EVEN BIOAEROSOLS COVER A BROAD RANGE



PHYSICS OF TRANSPORT, GROWTH

- Size (aerodynamic size?, shape is also important)
- Concentration (number based metric, $\# / \text{cm}^3$)
- Morphology – shape, agglomerate state
- Biological particles generally carry a surface charge, very important
- Other surface characteristics – hydrophobic, hydrophilic, etc.



EXPOSURE CHAMBER



AEROSOL GENERATOR

AEROSOL MEASUREMENT



OTHER IMP. CONSIDERATIONS

- Transport Modeling
- Aerosol dynamics (size change)
- Respiratory Deposition

Aerosol Generation

- Aerosolize suspension in liquid without changing viability (considerations: fragmentation, charge distribution, etc)
- Role of other constituents – surfactants, ionic salts, etc.
- Resultant throughput and size distribution

ATOMIZATION




- Uses a high velocity air jet to create a pressure difference that pulls liquid w/suspension from a reservoir
- Air jet breaks up liquid stream into droplets
- Impaction results in size selectivity – only smaller droplets are aerosolized. Larger droplets drop back into reservoir
- Possible disadvantage: damaging of organism on impaction. Re-design will result in gentler impaction, and non-reuse of larger droplets generated

NEBULIZERS


- Similar principle as atomizer – gentler aerosolization
- Can have higher throughputs, lesser control on size
- Designs (such as BANG, etc to minimize foaming)



SEM images of MS2 bacteriophage (polio virus surrogates) aerosolized with a Collison Nebulizer



(a) Virion particles in suspended solution



(b) Virion particles after nebulizer

Nebulized Particles Tend to Clump Together !

Hogan, Lee, Biswas (2003) In review, *Aerosol Sci. Technol.*


Though the viral particles are 25 nm in diameter, we measure a broader size distribution

Smaller particles are due to fragments of the particles – either due to nebulization or due to breakup as a result of surface charges

Larger particles due to tendency to aggregate

Our systems have allowed us to conduct fundamental studies on charge distributions on these viral particles – NEVER DONE BEFORE.



Why are such fundamental studies important? – To better design control methodologies, scale up, ensure byproducts are not formed




REAL TIME MEASUREMENT OF VIRAL PARTICLE SIZE DISTRIBUTIONS.

ELECTROSPRAYS

- Use an applied electrical field to break up droplets
- Ability to manipulate size in very narrow ranges and hence state of the resultant particles (overcome van der Waals forces which tend to hold particles together causing clumping)
- Can eliminate role of other additives in the solutions



POST GENERATION CONDITIONING



Summary of “Aerosol Generation”

- Nebulizers may not be able to produce single units of the bioagent; even if starting material is “ultra-pure”. Atomizers maybe better in size control, but may alter viability
- Post conditioning is a viable means of obtaining a well characterized aerosol, and maybe essential
- Electrosprays one of the best ways of generating particles in a large range of sizes

Aerosol Measurement

- Real time instruments that measure size distributions
- Most reliable are the electrical mobility techniques – especially for particles in the submicrometer and nanometer size ranges
- Optical instruments based on light scattering are also feasible – need some prior knowledge of organisms
- As real time instruments are well developed, routine monitoring in chambers using probes are advisable

GOVERNING EQUATIONS OF PARTICLE MOTION

- Equations for transport
- External force fields – gravity, electrical forces
- Other modes of transport: diffusion
- Important non-dimensional parameters

$$m \frac{d\bar{v}_p}{dt} = \frac{3\pi\mu d_p}{C} (\bar{u} - \bar{v}_p) + \bar{F}_{ext}$$


$$F_{electrical} = q E$$

$$D = \frac{kT}{3\pi\mu d_p} ; \text{Diffusion Coefficient}$$

$$St = \frac{\rho_p C d_p^2 U_0}{18\mu L}$$

DIFFERENTIAL MOBILITY ANALYZERS

- Differential Mobility Analyzer – known charge on particles that are then classified in an electrical field
- Based on tuned electrical field (applied voltage) – a very narrow size of the particles can be selected
- The monodisperse stream can be used for further studies, or can be sent to a Condensation Particle Counter for counting
- Sequential change in voltage allows mapping the entire size distribution – Instrument called the Scanning Mobility Particle Sizer (SMPS)



Overall system layout (SMPS)




Washington University in St. Louis

AA QRL - Aerosol and Air Quality Research Laboratory


NANO - DMA

- Need to measure and classify real small particles – nanometer sized, such as viruses, protein fragments, etc
- Diffusion losses preclude effective measurement
- Modify DMA to minimize flight time and losses of nanoparticles – the NANO- DMA.



INSTRUMENTATION FOR QUICK DETECTION OF VIRUSES & PROTEIN FRAGMENTS

MS- 2 BACTERIOPHAGE (POLIO VIRUS SURROGATE)




Routinely used atomizers
And DMA

Improved accuracy,
and reduced time,
Using Electrospray
System
(Kulkarni, Chen, Biswas, 2003)

System can also measure
Protein fragments - selectivity

**APPLY TO BIOAGENT
DETECTION SYSTEMS**

OPTICAL MEASUREMENT SYSTEMS




Schematic of Experimental Setup



Washington University in St. Louis

AA QRL - Aerosol and Air Quality Research Laboratory


Rayleigh and Mie Scattering Theories




$$I(\theta) = \frac{I_0 \lambda^2 (C_{vv} + C_{hh})}{8 \pi^2 R^2}$$

Scattered Intensity is a function of particle size, refractive index, angle of detection


MULTIANGLE LIGH SCATTERING MEASUREMENTS





Washington University in St. Louis

AA QRL - Aerosol and Air Quality Research Laboratory




Washington University in St. Louis

AA QRL - Aerosol and Air Quality Research Laboratory


Detection of Micro-Organisms

- There is a renewed interest in quick and robust methods of detection of micro-organisms
- Laser based optical scattering systems offer a potential method of identifying larger bacteria – depending on shape, size and viability, they have unique light scattering patterns as a function of angle
- Can also modify to measure uv-vis absorption, or fluorescent signal to improve selectivity of detection.

ANGULAR DISTRIBUTION OF SCATTERED INTENSITY FOR ROD SHAPED MICRO-ORGANISMS



ANGULAR DISTRIBUTION OF SCATTERED INTENSITY – OTHER MICRO-ORGANISMS



**DETECTION OF BIOLOGICAL
PARTICLES**




Washington University in St. Louis

AA QRL - Aerosol and Air Quality Research Laboratory

AGIs – How effective are they?

$$St = \frac{\rho_p C d_p^2 U_0}{18 \mu L}$$




- AGI – Impinger designed for collection based on impaction of bioaerosols onto liquid surface
- Impaction efficiency increases with Stokes Number, St; which is proportional to square of diameter
- Bacteria which is in 1 to 10 μm range is 100 times larger than viruses. Hence St is 10,000 times smaller. Collection Efficiency is not going to be 100%
- Increasing velocity may help – but may alter viability of organism during collection
- Need to design more effective collection systems for viruses and fragments

Summary of Aerosol Measurement Devices

- Total particle counts – best instrument is the Condensation Particle Counter (CPC)
- Particles less than 100 nm – use the NanoDMA – CPC system
- Particles between 100 and 1000 nm – use the DMA (SMPS system)
- Particles greater than 1000 nm (1 μm) – use Optical Particle Counters (OPC) or Aerodynamic Particle Sizers
- Characterize systems such as AGI's for the system in which they are to be used

Exposure chambers



- Is measuring the viable organisms by AGI's at the outlet good enough to calculate dose??

- Due to complex flow patterns there is a spatial variation of aerosol concentration, and hence may have localized regions – especially near animal intake ports

- Have the tools to characterize this – both experimentally and theoretically!



Washington University in St. Louis

AA QRL - Aerosol and Air Quality Research Laboratory

Flow modeling, distribution in chamber


- Governing equation for particle transport

$$\frac{\partial n}{\partial t} = \underbrace{D \nabla^2 n}_{\text{Diffusion}} - \underbrace{\bar{\mathbf{u}} \cdot \nabla n}_{\text{Convection}} + \underbrace{\frac{D}{kT} \nabla \cdot (\bar{\mathbf{F}}_{\text{EXT}} n)}_{\text{External Forces}}$$

- Need flow field – can solve Navier Stokes equations, or use CFD codes
- Use velocity field, \mathbf{u} , to solve particle size distribution, n , variation as a function of space and time
- Particle spatial distribution will be a function of size, geometry of chamber

Respiratory deposition


- Using spatial and time variation of size distribution function in respiratory chamber, can calculate size dependent deposition in respiratory system
- Importance of particle morphology (aggregates) and charge distribution very critical for accurate estimation
- <http://www.aerosols.wustl.edu/aaqr/courses/cycopcr esp/>




BIOAGENT INACTIVATION STUDIES: DEVELOPMENTAL PROJECT

- **Capture and inactivation of bioagents, such as viruses and bacteria, from air streams important in many applications**
 - **Applications in aircraft cabin air filtration, indoor building ventilation systems, chem-bio agent inactivation (counter measures to bioterrorism)**
- **Several agents of interest – SARS virus, Small pox virus, spores, anthrax**
- **Several technologies in various stages of development (Alving, 2002). Most are based on filtration systems which are plagued with high pressure drops, maintenance and operational problems**


NOVEL CORONA-SOFT XRAY SYSTEM FOR HIGH EFFICIENCY PARTICLE CAPTURE





- High voltage generates a corona, that produces an ion rich environment. Due to irradiation from Soft x-ray unit, there is a cascading effect that results in much higher ion concentrations
- Difficult to charge particles (ultrafine sizes) - are readily charged, and trapped in electrical field
- Due to high ion concentrations, and oxidizing environment, organic species are readily converted to carbon dioxide.
- Can further promote this with an nanostructured catalyst coating
- Very compact unit
- Can be readily mounted on existing duct work
- More than "4 to 5 logs" removal demonstrated for surrogates such as polio virus


Currently funded by the NIH Bio-Defense Project



Washington University in St. Louis


AA QRL - Aerosol and Air Quality Research Laboratory

Figure 2



EXPERIMENTAL SYSTEM FOR AEROSOLIZING VIRUSES FOR CAPTURE TESTS IN CORONA-SOFT X_RAY UNIT


Capture efficiency at different particle diameters from 10 – 225nm without neutralization with a Po-210 bipolar charger(Q=10lpm).



Washington University in St. Louis

AA QRL - Aerosol and Air Quality Research Laboratory


Capture efficiency at different particle diameters from 10 – 225nm with a Po-210 bipolar charger($Q=10$ hpm).




Washington University in St. Louis

AA QRL - Aerosol and Air Quality Research Laboratory

LOG CAPTURE OF VIRAL PARTICLES AS A FUNCTION OF APPLIED VOLTAGE





Aerosol aspects of aerobiology quite complex; but aerosol science and engineering is sufficiently advanced to address and resolve these issues.

Collaboration with main stream aerosol scientists strongly recommended

- **Thank you for your attention**
- **Questions?**

Aerosol and Air Quality Research Laboratory

Washington University in St. Louis

www.aerosols.wustl.edu/aaql

Particle Formation and Capture

SYNTHESIS & PROCESSING

- Aerosol Flame Reactors
- Aerosol Furnace Reactors
- Pristine and Doped TiO₂
- Nanostructured Coatings
- Magnetic Oxides

CONTROL & CAPTURE

- Soft X-Ray/Unipolar Coronas for Particle Charging & Capture
- Capture in Magnetic Fields
- Photocatalysis using TiO₂ – inactivation of bioaerosols and organic species
- Oxygen Enriched Coal Combustion/Hg Capture

Ambient Aerosols and Air Quality

- Diesel Engine Exhausts & Children's Health
- Morphology of Ambient PM
- Health Effects of Nanoparticles

Colloids/Bio-colloids

- Detection of Micro-organisms
- Electrostatic Filtration
- Inactivation of bio-agents

Sponsors

**USEPA, NSF, NIEHS, NIH,
Ameren, Hamamatsu**

Graduate Students

Kuk Cho; Shaohua Hu; Pramod Kulkarni; Prakash Kumar; Rafael McDonald; Ayano Niwa

Marina Smallwood; Achariya Suriyawong

Undergraduate Students Post Doctoral Fellow

David Weingelst; Damien Brewer; Dr. Miles

Dr. Pratim Biswas

Stifel and Quinette Jens Professor
Director, Environmental Engr. Sci.

Campus Box 1180, WUSTL
St. Louis, MO 63130

Email: pbiswas@seas.wustl.edu

Tel: 314-935-5482

Aerosol and Air Quality Research Laboratory

AERODYNAMIC PARTICLE SIZER

- Particles are accelerated in a flow field
- Velocity is a function of the Stokes Number
- Use two laser system to time movement along a fixed distance, and then back calculate particle size
- Detection method good for large particles, as Stokes number goes as particle size². Can measure > 1 micrometer very effectively



Washington University in St. Louis

AA QRL - Aerosol and Air Quality Research Laboratory