System Description CARIM

# Clinical Agent Repository Inventory Management Computer Database System

# **CARIM 2.2**



**System Description** 

NIAID/DMID Office

Washington, DC April 2007

# **Table of Contents**

Table of Figures	6
Introduction	8
Project Description	9
Background	9
Solution	9
System Overview	10
Functionality	
Security	10
Navigation	10
Components	10
Common	
Reports	11
Forms	11
Naming Conventions	11
Functional Process Description	12
Login	12
Receiving Module	
Receiving Module Dashboard	
Receiving Normal Inventory of Clinical Agents	
Normal Receiving Report	14
Normal Receiving Form	14
Return Drug Receipts	16
Return Drug Report	16
Return Drug Form	17
Shipping Module	18
Shipping Module Dashboard	
DMID IND For Human Use Shipments	
DMID IND For Human Use Shipment Report	19
DMID IND For Human Use Shipment Form	
Non-DMID IND For Human Use Shipments	
Non DMID For Human Use Shipment Report	21

Non DMID IND For Human Use Shipment Form	21
DMID IND Not For Human Use Shipments	22
DMID IND Not For Human Use Shipment Report	23
DMID IND Not For Human Use Shipment Form	23
Return Shipments	24
Return Shipment Report	24
Return Shipment Form	25
HSROAD Orders	26
HSROAD Orders Report	27
HSROAD Import Form	27
Inventory Module	29
Inventory Adjustment	29
Internal Transfer Drug Inventory	30
Reports Module	32
Receiving Reports	
Receiving Normal Report	
Receiving Return Report	
Shipping Reports	
DMID For Human Use Shipment Report	
NON-DMID For Human Use Shipment Report	
DMID Not For Human Use Shipment Report	
Return Shipment Report	
Inventory Reports	
Inventory Normal Report	
Inventory Return Used Report	
Inventory Return Unused Report	40
Inventory Quarantine Report	40
Inventory Summary Report	41
Return Report	
Inventory Adjustments Report	
Inventory Transfer Report	43
Product Reports	44
Transactions Report - Normal Inventory	44
Transaction Repot - Return used	45
Transaction Report - Return Unused	46
Transaction Report - Quarantined	16

Product Protocols Report	47
Summary Reports	48
Receiving Normal Summary Report	48
Receiving Return Summary Report	48
Shipping Summary Report	49
Alerts	49
Alert Log	50
System Administration Module	51
Carrier	51
CARIM Carriers Report	51
Add a New Carrier Form	52
Manufacturer	53
CARIM Manufacturers Report	53
Add a Manufacturer Form	53
Products	54
CARIM Products List Report	54
Create a New Product Form	55
Users	56
CARIM Users Report	56
Create a new User Form	56
Utilities	57
Appendix – Functional Process Flowcharts	59

# **Table of Figures**

Figure 1 CARIM Login Page	
Figure 2 Receiving Module Dashboard	
Figure 3 Normal Receiving Report	
Figure 4 Normal Receiving Form	
Figure 5 Return Drug Report	
Figure 6 Return Drug Form	
Figure 7 Shipping Module Dashboard	
Figure 8 DMID IND For Human Use Shipment Report	
Figure 9 DMID For Human Use Shipment Form	
Figure 10 Non-DMID IND For Human Use Shipment	
Figure 11 Non DMID IND For Human Use Shipment Form	
Figure 12 DMID Not For Human Use Shipments	
Figure 13 DMID IND Not For Human Use Shipment Form	
Figure 14 Return Shipment Report	
Figure 15 Return Drug Shipment Form	
Figure 16 HSROAD Order Report	27
Figure 17 HSROAD Order Import Form	
Figure 18 Inventory Module	
Figure 19 Inventory Adjustment FormFigure 20 Inventory Transfer Form	
Figure 21 Reports Module	
Figure 22 Receiving Normal Report	
Figure 24 DMID For Human Use Shipment Report	
· · · ·	
Figure 25 Non DMID For Human Use Report	
Figure 27 Return Shipment Report	
Figure 28 Normal Inventory Report	
Figure 29 Inventory Return Used Report	
Figure 30 Inventory Return Unused Report	
Figure 31 Inventory Quarantine Report	
Figure 32 Inventory Summary Report	
Figure 33 Return Drug Inventory Report	
Figure 34 Inventory Adjustment Report	
Figure 35 Inventory Transfer Report	
Figure 36 Product Transactions Report - Normal Inventory	
Figure 37 Product Transactions Report - Return Used Inventory	
Figure 38 Product Transaction Report - Return Unused Inventory	
Figure 39 Product Transaction Report - Quarantine Inventory	
Figure 40 Product Protocol Report	
Figure 41 Receiving Normal Summary Report	
Figure 42 Receiving Return Summary Report	
Figure 43 Shipping Summary Report	
Figure 44 Alert Log	
Figure 45 CARIM Carriers Report	
Figure 46 Add a Carrier Form	
Figure 47 CARIM Manufacturer's Report	
Figure 48 Add a Manufacturer Form	
Figure 49 CARIM Products List Report	
Figure 50 Create a New Product Form	
Figure 51 CARIM User Report	
g 5 . 5 . 4	50

# **CARIM**

Figure 52 Add a New User Form	57
Figure 53 CARIM Utilities	58

System Description **CARIM** 

#### Introduction

CARIM is an inventory management computer database system developed to effectively manage and control all investigational drugs at the Division of Microbiology and Infectious **Diseases Clinical Agents Repository (DMID-CAR).** 

The Clinical Agents Repository Inventory Management (CARIM) system has been designed and developed by Fisher Bioservices and Cascades Technologies Incorporated (CTI). CARIM is intended to manage all investigational drugs at the Division of Microbiology and Infectious Diseases Clinical Agents Repository (DMID-CAR).

There are three components that must be considered by users who want to access the WWW: hardware, software, and network connectivity. All three must be present and established in order to access the web effectively. CARIM is based upon World Wide Web (WWW) technologies. As such, it consists of a single database accessed by users across the Internet using a web browser, such as Netscape's Navigator or Microsoft's Internet Explorer. This approach provides the advantages of a centrally supported system as well as increased accessibility through the Internet. Using a web browser, users access the CARIM server and log into the application. After completing the login process, the user is presented with welcoming information and a menu of available options. Information related to CARIM data can then be entered, retrieved, or modified by the user just as if the data were located on his or her personal computer.

System Description CARIM

# **Project Description**

# **Background**

The CARIM computer database system has been designed to automate the receiving, shipping, inventory control, and activity report processes of the DMID Clinical Agents Repository (DMID-CAR).

#### Solution

CARIM has been developed with Internet usability in mind. The Web server software for this project is Microsoft's Internet Information Services 6.0 (IIS). The ColdFusion Enterprise MX application server is used as the application server to connect the Web application to the database and provide dynamic content for the browser based CARIM application. All this architecture is based on the Microsoft's Windows 2003 server software. Fisher BioServices Inc, is providing the necessary server hardware and software licenses required for the project and it is targeted for users with version 5.0 or above of the Microsoft's Internet Explorer or Netscape Navigator browsers. The screens have been developed for a minimum screen resolution of 800 x 600 pixels. This will ensure that users with a wide array of monitor sizes and resolutions can see the system screens without scrolling sideways. The database has been developed following best practices for relational database design, data normalization, and development. This database is deployed on the Microsoft's SQL Server 2000 database server.

System Description **CARIM** 

# **System Overview**

The DMID-CAR utilizes the CARIM application as the tool for receipt, processing, tracking and fulfillment of the products orders used in clinical trials. This section provides a brief overview of system features and functionality available in the CARIM application. Topics covered in this section include application functional overview, application security, navigation, and components.

# **Functionality**

- Receiving clinical agents from manufacturers and clinical sites
- Shipping clinical agents to manufacturers and clinical sites
- Adjustments and transfers of internal inventory
- System reports
- System administration

# Security

- 3 permission levels
- Log in protection
- Sessions are tracked

# **Navigation**

- Main / application menu on top
- Page / module menu on left side

#### Components

Pages have specific components visible by the type of page. Most components are application wide and usually visible on all pages. Reports have specific components only available to reports. Forms have specific components only available to forms.

#### Common

- Page Title
- Introduction
- Report or form
- Header
- Footer
- Main menu
- Page menu

# Reports

CARIM reports all contain page titles, report search boxes, tabular data, edit and delete controls, legend, pagination, user messaging, and an add record link

#### **Forms**

CARIM forms all contain page title, form controls to add or edit data, user messaging, and confirmation pages

# **Naming Conventions**

- Receipts (YYR-MMDDXX or YYRS-MMDDXX)
- Shipments (YYS-MMDDXX)
- HSROAD Orders (YYO-MMDDXX)
- Adjustments (YYA-MMDDXX)
- Transfers (YYT-MMDDXX)

The system features and its related functionality make this application an effective solution for automating the managing the Repository transaction process and creates a bridge for supporting the ordering products via the Order It Module in HSROAD. The following sections will provide a description of the functional processes of the application.

System Description CARIM

# **Functional Process Description**

The CARIM System is used by the DMID - CAR to conduct all repository transactions. These transactions include the receipt and shipment of investigational agents which are to be used in Clinical trials. The following sections will provide an overview of the features and functionality available within each module of the application which automates these processes.

# Login

When the user accesses CARIM from a web browser, they are presented with a log in screen. If an account has been pre-established, the user will use that User ID and password to gain access to the system. If no account has been provided, a user will not be able to access the application.

There is access privileges associated with the application. The application only displays the modules which the user has access to. If a particular module is inaccessible to the user, it won't be displayed. The modules that accessible to each user have been predefined by the system administrator and are linked to the user's account.



Figure 1 CARIM Login Page

# **Receiving Module**

The Receiving Module is used to create, modify or delete all information pertaining to all received drug shipments. It is the process used to log products received by the Clinical Agent Repository into the CARIM system. The receiving module consists of two sub-modules and a dashboard.

Adding receipt records to CARIM follows a two step workflow. First, a new receipt is created by the receiving inspector. The receiving inspector is the user who first creates the receipt record in the system. The same or a different user, also defined as a receipt inspector may edit the record as often as required. Once created (and perhaps edited), another different user is required to confirm the receipt. Once confirmed the normal balance of inventory is affected, and the record is read-only in the system.

#### **Receiving Module Dashboard**

The receiving module dashboard gives the user a single page view of all activity pending in the receiving module. From here all normal and return receipts are visible. Also, the user has links to view all receipts by type (normal or return) or to add a new receipt by type (normal or return).



Figure 2 Receiving Module Dashboard

#### **Receiving Normal Inventory of Clinical Agents**


This area is designed to view, create, modify, or delete normal drug receipts. Typically, normal drug receipts are defined as products received from the manufacturer and shipped to the Clinical Agent Repository, then entered into the CARIM application. A transaction within this sub-module will affect the Normal, Rejected and Quarantined inventories. Once receipts are confirmed, the inventory of that product increases.

These types of receipts are designated in the system by a series of numbers and a letter "R" to identify it as a receipt. A typical drug receipt number may be 02R-021201. These numbers equate to transactions within the inventory module. See the Naming Convention section for more information.

#### **Normal Receiving Report**

The normal drug receipt information report displays a list of all clinical agents received in CARIM in the past 3 years adding to the normal inventory balance. Completed transactions will have a confirm date visible; uncompleted transactions will not. From this report the user may add a new drug to the normal inventory, or edit or delete an existing but yet to be completed transaction. Completed transactions are read only. Users may also search all receipt transactions using the search box or the record set pagination controls.

Data displayed in this report includes the Receipt Number, Product number, Product Name (ID lot number – label name), the receipt inspector name, the receipt confirmation name, and the date the record was confirmed. An action column is displayed providing icons to edit or delete unconfirmed normal receipt records.




**Figure 3 Normal Receiving Report** 

#### **Normal Receiving Form**

The appropriate form to add, edit or delete a normal receipt is accessible by clicking the *Add normal drug receipt* link, clicking the edit icon, or clicking the delete icon next to the desired record.

This form displays the form controls required and optional defined for normal drug receipts. Required fields are denoted by the red line to the right of the form field label.



**Figure 4 Normal Receiving Form** 

#### **Return Drug Receipts**

This area is designed to view, create, modify, or delete return drug receipts. Typically, return drug receipts were previously shipped to Protocol Investigators at the clinical site and then returned to the Clinical Agent Repository. A transaction within this sub-module will affect the return used and return unused inventories. Once receipts are confirmed, the inventory of that product increases.

These types of receipts are identified in the system by a series of numbers and the letters "RR" to identify it as a return receipt. A typical return drug receipt number may be 02RR-031801. These associate to transactions within the inventory module. A transaction within this submodule will affect the Return Unused and Return Used repositories.

#### **Return Drug Report**

The return drug receipt information report displays a list of all clinical agents returned to the repository and logged CARIM in the past 3 years. Completed return transactions will have a confirm date visible; uncompleted transactions will not. From this report the user may add a new drug to the return used or unused inventories, or edit or delete an existing but yet to be completed transaction. Completed transactions are read only. Users may also search all receipt transactions using the search box or the record set pagination controls.

Data displayed in this report includes the Receipt Number, Product Name (ID - lot number – label name), the receipt inspector name, the receipt confirmation name, and the date the record was confirmed. An action column is displayed providing icons to edit or delete unconfirmed return receipt records.



Figure 5 Return Drug Report

#### **Return Drug Form**

The appropriate form to add, edit or delete a return receipt is accessible by clicking the *Add* return drug receipt link, clicking the edit icon, or clicking the delete icon next to the desired record.

This form displays the form controls required and optional defined for adding return drug receipts into CARIM. Required fields are denoted by the red line to the right of the form field label.




Figure 6 Return Drug Form

System Description CARIM

# **Shipping Module**

The Shipping Module is used to create, modify or delete all information pertaining to all product shipments. It is the process used to log the shipped products into the CARIM system.

#### **Shipping Module Dashboard**

The shipping module dashboard gives the user a single page view of all pending activity in the shipping module. From here all DMID IND for human use, Non-DMID for human use, DMID not for human use, return shipments, and HSROAD orders are visible. Also, the user has links to view all shipments by type or to add a new receipt by type.




Figure 7 Shipping Module Dashboard

#### **DMID IND For Human Use Shipments**

This area is designed to create, modify, or delete a shipment of DMID IND clinical agents intended for human use. These types of shipments are identified in the system by a series of numbers and the letter "S" to identify it as a shipment, a typical drug shipment number may be 02S-031801. For more information refer to the previous section - Naming Conventions. These in turn become transactions within the inventory module affecting only the Normal repository.

#### **DMID IND For Human Use Shipment Report**

The DMID IND for Human Use Shipment report displays a list of all For Human Use shipments made from the Repository to clinical sites in the past 3 years subtracting from the normal inventory balance. Completed transactions will have a confirm name visible; uncompleted transactions will not. From this report the user may add a new shipment, edit or delete an existing but yet to be completed transaction. Completed transactions are read only. Users may also search all receipt transactions using the search box or the record set pagination controls.

Data displayed in this report includes the Shipping Number, Shipping Date, the Shipping inspector name, the Shipping confirmation name, and the name of recipient Clinical Site. An action column is displayed providing icons to edit or delete unconfirmed DMID IND For Human Use records.




Figure 8 DMID IND For Human Use Shipment Report

#### **DMID IND For Human Use Shipment Form**

The appropriate form to add, edit or delete a shipment of this type is accessible by clicking the *Add a DMID IND for Human Use Shipment* link, clicking the edit icon, or clicking the delete icon next to the desired record.

This form displays the form controls required and optional as is defined for DMID IND For Human Use shipments. Required fields are denoted by the red line to the right of the form field label.




Figure 9 DMID For Human Use Shipment Form

#### **Non-DMID IND For Human Use Shipments**

This area is designed to create, modify, or delete a shipment of Non-DMID IND clinical agents intended for human use. These types of shipments are identified in the system by a series of numbers and the letter "S" to identify it as a shipment. A typical drug shipment number may be 02S-031801. Refer to the previous section - Naming Conventions. These in turn become transactions within the inventory module affecting only the Normal repository.

#### Non DMID For Human Use Shipment Report

The Non DMID IND for Human Use Shipment report displays a list of all shipments made from the Repository to clinical sites in the past 3 years subtracting from the normal inventory balance. Completed transactions will have a confirm name visible; uncompleted transactions will not. From this report the user may add a new shipment, edit or delete an existing but yet to be completed transaction. Completed transactions are read only. Users may also search all receipt transactions using the search box or the record set pagination controls.

Data displayed in this report includes the Shipping Number, Shipping Date, the Shipping inspector name, the Shipping confirmation name, and the Carrier Name. An action column is displayed providing icons to edit or delete unconfirmed Non DMID IND For Human Use records




Figure 10 Non-DMID IND For Human Use Shipment

#### Non DMID IND For Human Use Shipment Form

The appropriate form to add, edit or delete a shipment of this type is accessible by clicking the *Add a Non-DMID IND for Human Use Shipment* link, clicking the edit icon, or clicking the delete icon next to the desired record.

This form displays the form controls required and optional as is defined for Non DMID IND For Human Use shipments. Required fields are denoted by the red line to the right of the form field label.




Figure 11 Non DMID IND For Human Use Shipment Form

#### **DMID IND Not For Human Use Shipments**

This area is designed to create, modify, or delete a shipment of DMID IND clinical agents not intended for human use. These types of shipments are identified in the system by a series of numbers and the letter "S" to identify it as a shipment. A typical drug shipment number may be 02S-031801. For more information refer to the previous section - Naming Conventions. These in turn become transactions within the inventory module affecting only the Normal repository.

#### **DMID IND Not For Human Use Shipment Report**

The DMID IND Not for Human Use Shipment report displays a list of all Not For Human Use shipments made from the Repository to clinical sites in the past 3 years subtracting from the normal inventory balance. Completed transactions will have a confirm name visible; uncompleted transactions will not. From this report the user may add a new shipment, edit or delete an existing but yet to be completed transaction. Completed transactions are read only. Users may also search all receipt transactions using the search box or the record set pagination controls.

Data displayed in this report includes the Shipping Number, Shipping Date, the Shipping inspector name, the Shipping confirmation name, and the Carrier name. An action column is displayed providing icons to edit or delete unconfirmed DMID IND Not For Human Use records.




Figure 12 DMID Not For Human Use Shipments

#### **DMID IND Not For Human Use Shipment Form**

The appropriate form to add, edit or delete a shipment of this type is accessible by clicking the *Add a DMID IND Not for Human Use Shipment* link, clicking the edit icon, or clicking the delete icon next to the desired record.

This form displays the form controls required and optional as is defined for DMID IND Not For Human Use shipments. Required fields are denoted by the red line to the right of the form field label.




Figure 13 DMID IND Not For Human Use Shipment Form

#### **Return Shipments**

This area is designed to create, modify, or delete a shipment of clinical agents to be returned to their manufacturer. These types of shipments are identified in the system by a series of numbers and the letters "RS" to identify it as a return shipment. A typical return drug shipment number may be 02RS-031801. For more information refer to the previous section - Naming Conventions. A transaction within this sub-module will affect the Return Unused and Return Used inventory.

#### **Return Shipment Report**

The Return Shipment report displays a list of all Return drug shipments returned to the manufacturer in the past 3 years affecting the Return Unused and Return Used inventory balance. Completed transactions will have a confirm name visible; uncompleted transactions will not. From this report the user may add a new shipment, edit or delete an existing but yet to

be completed transaction. Completed transactions are read only. Users may also search all receipt transactions using the search box or the record set pagination controls.

Data displayed in this report includes the Shipping Number, Shipping Date, the Shipping inspector name, the Shipping confirmation name, and the Carrier name. An action column is displayed providing icons to edit or delete unconfirmed Return Drug shipments.




**Figure 14 Return Shipment Report** 

#### **Return Shipment Form**

The appropriate form to add, edit or delete a shipment of this type is accessible by clicking the *Add a Return Drug shipment link*, clicking the edit icon, or clicking the delete icon next to the desired record.

This form displays the form controls required and optional as is defined for Return Drug shipments. Required fields are denoted by the red line to the right of the form field label.



**Figure 15 Return Drug Shipment Form** 

#### **HSROAD Orders**

This area is designed to as a starting point to fulfill clinical agent orders placed in HSROAD. It is here that orders placed in the HSROAD application are imported into the CARIM application. Once a HSROAD order has been imported, it will be available through the appropriate option based on the transaction type. At this time, the only supported order type that may be placed/fulfilled from HSROAD are DMID IND clinical agents intended for human use trials.

For more information on ordering agents for clinical trials, please refer to the HSROAD system description section on the Order It Module, or the Order It Module System description.

#### **HSROAD Orders Report**

The HSROAD Orders report displays a list of all orders pending import into the system for fulfillment. Transactions of this type affect the normal inventory balance. The only action available in this report is the Import option, once imported these orders are removed from this report and are listed with other transactions in the DMID For Human Report. They will continue to be added by their HSROAD order number. Please refer to the previous section - Naming Conventions for more information. Users may also search all receipt transactions using the search box or the record set pagination controls.

Data displayed in this report includes the Order Number, Protocol Number, Product (ID - Lot-Label); Order Creator, and the Date Ordered. An action column is displayed providing an icon to import the HSROAD order.



Figure 16 HSROAD Order Report

#### **HSROAD Import Form**

Clicking the *Import* icon next to the desired record displays the Import HSROAD Order Form. This form displays the order information from the HSROAD order form, related comments, Recipient contact information and required form controls as is defined for HSROAD Order shipments. Required fields are denoted by the red line to the right of the form field label.





Figure 17 HSROAD Order Import Form

System Description CARIM

# **Inventory Module**

The Inventory module is used to adjust or transfer the balance of units available from each type of inventory tracked by CARIM. This module can be entered from the Inventory tab on the CARIM main menu (for users that have access to the module). Once the tab is clicked, the user will be presented with a page/module-level menu to each of the individual sub-modules.



**Figure 18 Inventory Module** 

#### **Inventory Adjustment**

This sub-module was designed to make adjustments to the current inventory tables if an error has occurred and needs to be modified from the original entry. The inventory adjustment module was created in order to avoid changes being made directly on the database. Once a record has been confirmed within the system it cannot be modified. The inventory adjustment form is then useful in order to make the necessary adjustment.


The Inventory adjustment form displays the form controls both required and optional as is defined for making this type of adjustment. Required fields are denoted by the red line to the right of the form field label.



**Figure 19 Inventory Adjustment Form** 

#### **Internal Transfer Drug Inventory**

This sub-module was designed to transfer drug units from one inventory table to another within the database. Sometimes a product has been logged into the wrong repository or its classification needs to be changed. Instead of an adjustment, a simple transfer is in order. Via the transfer form it is as simple as designating the product to transfer, the repository it is being transferred from, the repository it is going to, and the units to be transferred.



**Figure 20 Inventory Transfer Form** 

System Description CARIM

# **Reports Module**

The Report module is used to display reports of the information from the CARIM database, such as receivables, shipments, inventory, etc. This module can be entered from the CARIM main menu for users that have access to it. The user will be presented with a links menu to each of these individual reports group by categories. Those categories are:

- ✓ Receiving
- ✓ Shipping
- ✓ Inventory
- ✓ Product
- ✓ Summary




Figure 21 Reports Module

Each report has features such as searches, sorting, and printing. The search box allows for searches by a text string or by date range. Reports in tabular format can be sorted by clicking on column headers. Clicking on the Print Preview Button will take the user to a printable version of the report, which will not have any of the navigation bars and will have all records for the report on a single page. The user will have the ability to print their selected report via the browser window print option.


#### **Receiving Reports**

These reports provide historical information on all transactions related to Normal and Return receipts. The system maintains three years of historical data which can be searched by text, dates or product labels within the report. Detailed information on any receipt can be viewed by clicking the Receipt Number link.

#### **Receiving Normal Report**

The Receiving Normal report lists the normal receipts and their respective information. The report displays the following information fields per Receipt number:

- ✓ Receipt Number
- ✓ Confirmed Date
- ✓ Product (Product ID Lot Number Label Number)
- ✓ Volume per Container
- ✓ Received
- Quarantined
- √ Rejected



**Figure 22 Receiving Normal Report** 

#### **Receiving Return Report**

The Receiving Return report lists the returned receipts and their respective information. The report displays the following information fields per Receipt number:

- ✓ Receipt number
- ✓ IND Number/Protocol Number
- ✓ Protocol Investigator/Site

- ✓ Product ( Product ID Lot Number Label Number)
- ✓ Volume per Container
- ✓ Units Returned Used
- ✓ Units Returned Unused
- Units Quarantined



Figure 23 Receiving Return Report

#### **Shipping Reports**

These reports provide historical information on all shipments made by Product type. The system maintains three years of historical data which can be searched by text, dates or product labels within the report. Detailed information on any shipment can be viewed by clicking the Ship Number link.

#### **DMID For Human Use Shipment Report**

The DMID for Human Use Shipment report lists the DMID IND shipments for Human Use agents and their respective information. The report displays the following information fields per shipment number:

- ✓ Shipment Number
- ✓ Confirm Date
- ✓ IND Number/Protocol Number
- ✓ PI /Performance Site
- ✓ Units Shipped
- ✓ Product (Product ID Lot Number Label Number
- ✓ Unit Type



**Figure 24 DMID For Human Use Shipment Report** 

## **NON-DMID For Human Use Shipment Report**

The NON-DMID for Human Use Shipment report lists the NON-DMID shipments for Human Use agents and their respective information. The report displays the following information fields per shipment number:

- ✓ Shipment Number
- ✓ Confirm Date
- ✓ PI /Performance Site
- ✓ Product (Product ID Lot Number -Label Name)
- ✓ Units Shipped
- ✓ Unit Type



Figure 25 Non DMID For Human Use Report

# **DMID Not For Human Use Shipment Report**

The DMID Not for Human Use shipment report lists the DMID shipments Not For Human Use agents and their respective information. The report displays the following information fields per return shipment number:

- ✓ Shipment Number
- ✓ Confirm Date
- ✓ PI /Performance Site
- ✓ Product (Product ID Lot Number -Label Name)
- ✓ Units Shipped
- ✓ Container



Figure 26 DMID IND Not For Human Use Shipment Report

#### **Return Shipment Report**

The Return Shipment report lists the DMID returned shipments to manufacturers and their respective information. The report displays the following information fields per return shipment number:

- ✓ Shipment Number
- ✓ Confirm Date
- ✓ IND Number/Protocol Number
- ✓ Product (Product ID Lot Number Label Name)
- ✓ Manufacturer
- ✓ Units Returned Used
- ✓ Units Returned Unused
- ✓ Unit Type



**Figure 27 Return Shipment Report** 


#### **Inventory Reports**

These reports provide historical information on all Inventory adjustments and transfers made to all repository balance for the past three years. Reports are categorized by repository type for easy access. These reports can be searched by text and product label via the search option available in each report.

# **Inventory Normal Report**

The Inventory Normal Report lists the balances for all Normal product receipts in the repository. The report displays the following information fields:

- ✓ Product (Product ID Lot Number Label Name)
- ✓ Manufacturer
- ✓ Unit Type
- ✓ Current Balance
- ✓ Expiration Date



**Figure 28 Normal Inventory Report** 

# **Inventory Return Used Report**

The Inventory Return Used Report lists the balances for all returned used products in the repository. The report displays the following information fields:

- ✓ Product (Product ID Lot Number Label Name)
- ✓ Manufacturer
- ✓ Unit Type
- ✓ Current Balance
- ✓ Expiration Date




Figure 29 Inventory Return Used Report

# **Inventory Return Unused Report**

The Inventory Return Unused Report lists the balances for all returned unused products in the repository. The report displays the following information fields:

- ✓ Product (Product ID Lot Number Label Name)
- ✓ Manufacturer
- ✓ Unit Type
- ✓ Current Balance
- ✓ Expiration Date




**Figure 30 Inventory Return Unused Report** 

# **Inventory Quarantine Report**

The Inventory Quarantine Report lists the balances for all products in quarantine in the repository. The report displays the following information fields:

- ✓ Product (Product ID Lot Number Label Name)
- ✓ Manufacturer
- ✓ Unit Type
- ✓ Current Balance
- ✓ Expiration Date




**Figure 31 Inventory Quarantine Report** 

# **Inventory Summary Report**

The Inventory Summary Report lists the balances for all products in the repositories. The report displays the following information fields:

- ✓ Product (Product ID Lot Number Label Name)
- √ Normal Balance
- ✓ Return Used Balance
- ✓ Return Unused Balance
- ✓ Quarantine Balance
- ✓ Current Total Balance
- ✓ Unit Type
- ✓ Expiration Date



**Figure 32 Inventory Summary Report** 

#### **Return Report**

The Return Drug Inventory Report lists the balances for all returned and rejected products in the repository. The report displays the following information fields:

- Product (Product ID Lot Number Label Name)
- Manufacturer
- Returned Balance
- Rejected Balance
- **Expiration Date**



**Figure 33 Return Drug Inventory Report** 

# **Inventory Adjustments Report**

The Inventory Adjustments Report lists all adjustments made within the repository. The report displays the following information fields: This report can be searched by text, From and To dates and also Product Label.

- ✓ Adjustment Number
- ✓ Adjustment Date
- Product (Product ID Lot Number Label Name)
- Manufacturer
- **Expiration Date**
- ✓ Inventory
- Units Available
- Units Adjusted
- **New Balance**



**Figure 34 Inventory Adjustment Report** 

# **Inventory Transfer Report**

The Inventory Transfer Report lists all transfers made between the repositories. The report displays the following information fields:

- ✓ Adjustment Number
- √ Adjustment Date
- ✓ Product (
- ✓ Manufacturer
- ✓ Expiration Date
- √ From Inventory
- √ To Inventory
- ✓ Units Transferred
- ✓ Comments



**Figure 35 Inventory Transfer Report** 

#### **Product Reports**

The Product Reports provided historical information on all confirmed transactions by transaction type. Transaction reports are provided for the Normal Return Used, Return Unused and Quarantine Inventory transactions and also Product reports by Protocol and Expiration dates. Each report can be search by Text, From and To dates and Product label.

# **Transactions Report - Normal Inventory**

The Product Transactions Report - Normal Inventory lists all the confirmed transactions which affected the Normal Inventory in CARIM. The report displays the following information fields:

- ✓ Transaction Date
- ✓ Transaction Number
- ✓ Transaction Type
- ✓ Product (Product ID Lot Number Label Name)
- ✓ Units
- √ Balance




Figure 36 Product Transactions Report - Normal Inventory

## **Transaction Repot - Return used**

The Product Transactions Report - Return Used lists all the confirmed transactions which affected the Return Used Inventory in CARIM. The report displays the following information fields:

- ✓ Transaction Date
- ✓ Transaction Number
- ✓ Transaction Type
- ✓ Product (Product ID Lot Number Label Name)
- ✓ Units
- √ Balance



Figure 37 Product Transactions Report - Return Used Inventory

# **Transaction Report - Return Unused**

The Product Transactions Report - Return Unused lists all the confirmed transactions which affected the Return Unused Inventory in CARIM. The report displays the following information fields:

- Transaction Date
- **Transaction Number**
- Transaction Type
- Product (Product ID Lot Number Label Name)
- Units
- Balance



Figure 38 Product Transaction Report - Return Unused Inventory

## **Transaction Report - Quarantined**

The Product Transactions Report - Quarantined Inventory lists all the confirmed transactions which affected the Quarantine Inventory in CARIM. The report displays the following information fields:

- ✓ Transaction Date
- ✓ Transaction Number
- **Transaction Type**
- Product (Product ID Lot Number Label Name)
- Units
- Balance




Figure 39 Product Transaction Report - Quarantine Inventory

# **Product Protocols Report**

The Product Protocol Report lists all the confirmed transactions by Protocol number/Protocol short Name. The report displays the following information fields:

- ✓ Product (Product ID Lot Number Label Name)
- **Expiration Date**
- Protocol Number/Short Name



**Figure 40 Product Protocol Report** 


# **Summary Reports**

The Summary Reports provide quick reference to summary totals for Return, Return Used and Shipped transactions. Each report provides totals with no links for accessing detailed information. Transactions can be searched by date range.

# **Receiving Normal Summary Report**

The Receiving Normal Summary report lists the total of normal receipts. The report displays the following information fields:

- ✓ Receipts Total
- ✓ Total Units Received
- ✓ Total Units Quarantined
- ✓ Total Units Rejected



**Figure 41 Receiving Normal Summary Report** 

#### **Receiving Return Summary Report**

The Receiving Return Summary report lists the total of returned receipts. The report displays the following information fields:

- ✓ Receipts Total
- ✓ Total Units Received.
- ✓ Total Used
- ✓ Total Units Quarantined





Figure 42 Receiving Return Summary Report

## **Shipping Summary Report**

The Shipping Summary report lists totals of all shipments for investigational use, not for human use and return drug types. The report displays the following information fields:

- ✓ Shipment Type
- ✓ Total Shipments
- ✓ Total Sites
- ✓ Total Units Shipped



**Figure 43 Shipping Summary Report** 

#### **Alerts**

Alerts or Email notifications are generated when certain events occur within the systems. These events include advices sent to Order Creators and associated personnel for performance site when their placed orders have been confirmed and shipped to the performance site. Scheduled nightly task generate Product Expiry email alerts which advise system users of those products

which will expire in either 60 or 42 days. Archive email alerts are also included in the nightly task run and will archive all transactions which meet the following criterion:

Archive Products, Shipments and Receipts Criteria:

- ✓ Products with 0 balance for more than 180 days
- ✓ Shipments which have been confirmed for more than 1095 days.
- ✓ Receipts confirmed for more than 1095 days

# **Alert Log**

The alert log report displays a list of all email notifications sent from the CARIM application. The report displays the following information fields:

- ✓ Alert #
- ✓ Sent To
- ✓ Sent From
- ✓ Subject
- ✓ Attachment
- ✓ Date of Attempted Send

Clicking the subject of each report allows the user to view the details of the specific alert chosen. The details will display all information retained about the email notification, including a recipient list, subject, body, name of any attachments, format, and attempted send date.





Figure 44 Alert Log

System Description CARIM

# **System Administration Module**

The System Administration Module is used to create, edit, or delete the carriers, manufacturers, products, users, and look-up table values used throughout the system. This module can be entered from the CARIM main menu for users that have administrative access. The user will be presented with a links menu to each of the individual sub-modules for the system administration module.




#### Carrier

This sub-module is used to create, modify, or delete the carrier's information for the drugs received or shipped within the CARIM system. The carriers entered via this sub-module populate the drop-down lists within the application.

#### **CARIM Carriers Report**

The CARIM Carriers report displays a list of all Carriers currently entered into the system. Users may also search records using the search box or the record set pagination controls.

Data displayed in this report includes the Action, Carrier, Address; Contact Name, Phone number, Fax number, Email and Archived. An action column is displayed providing icons to edit or delete entered records.



**Figure 45 CARIM Carriers Report** 

#### Add a New Carrier Form

The appropriate form to add, edit or delete a Carrier record is accessible by clicking the *Add a New Carrier link*, clicking the edit icon, or clicking the delete icon next to the desired record.

This form displays the form controls required and optional as is defined for adding Carriers to the system. Required fields are denoted by the red line to the right of the form field label.



Figure 46 Add a Carrier Form

#### Manufacturer

This sub-module is used to create, modify, or delete the manufacturer's information for the drugs received or shipped within the CARIM system. The manufacturers entered via this sub-module populate the drop-down lists within the application.

# **CARIM Manufacturers Report**

The CARIM Manufacturers report displays a list of all Manufacturers entered into the system. Users may also search records using the search box or the record set pagination controls.

Data displayed in this report includes the Action, Manufacturer, Address; Contact Name, Phone number, Fax number, Email and Archived. An action column is displayed providing icons to edit or delete entered records.




Figure 47 CARIM Manufacturer's Report

#### Add a Manufacturer Form

The appropriate form to add, edit or delete a Manufacturer record is accessible by clicking the *Add a New Carrier link*, clicking the edit icon, or clicking the delete icon next to the desired record.

This form displays the form controls required and optional as is defined for adding a Manufacturer to the system. Required fields are denoted by the red line to the right of the form field label.



Figure 48 Add a Manufacturer Form


#### **Products**

This sub-module is used to create, modify, or delete all received or shipped product information within the CARIM system. Products are associated with the respective protocols and made accessible to HSROAD users via the Order It module with the Display in HSROAD options. Products with this option set to No will not be visible in the product list for a particular protocol. For more information on Ordering products through the Order it Module please refer to the either the Order it Module section of the HSROAD System Description or the Order it Module System Description. The products entered via this sub-module populate the drop-down lists within the application.

#### **CARIM Products List Report**

The CARIM Product List report displays a list of all products entered into the system both in use and archived. Users may also search records using the search box or the record set pagination controls.

Data displayed in this report includes the Action, Product (Product ID - Lot Number - Label Name), Container (Volume Weight/Unit Type); Manufacturer Name, Expiration Date, Archived, and Display in HSROAD. An action column is displayed providing icons to edit or delete entered records.



**Figure 49 CARIM Products List Report** 

#### **Create a New Product Form**

The appropriate form to add, edit or delete a Product record is accessible by clicking the *Add a New Product link*, clicking the edit icon, or clicking the delete icon next to the desired record.

This form displays the form controls required and optional as is defined for adding Products to the system. Required fields are denoted by the red line to the right of the form field label.




Figure 50 Create a New Product Form

System Description CARIM

#### Users

This sub-module is used to create, modify, or delete all pertinent information for the users with access to the CARIM system. Those users with access to this module have the ability to reset user passwords and send an alert to the affected user advising of the change. System permissions are set in this module and consist of three levels. These levels include:

- ✓ Read Only
- ✓ Read & Write
- ✓ Reports

## **CARIM Users Report**

The CARIM Users report displays a list of all users and user account information that have access to the system. Users may also search records using the search box or the record set pagination controls.

Data displayed in this report includes the Action, User ID, Name; E-mail, Last Login Date, Read Only, Read & Write and Reports. An action column is displayed providing icons to edit or delete existing records.




Figure 51 CARIM User Report

#### Create a new User Form

The appropriate form to add, edit or delete a user record is accessible by clicking the *Create a new user link*, clicking the edit icon, or clicking the delete icon next to the desired record.

This form displays the form controls required and optional as is defined for adding Products to the system. Required fields are denoted by the red line to the right of the form field label.




Figure 52 Add a New User Form

#### **Utilities**

The Utilities module is used to edit the Lookup tables for the application. This module can be entered from the System Administration Module for users with administrative access. Upon entering the Utilities module, the administrator will be presented with a table listing all the available for edit lookup tables in the system. The tables can be sorted by the table name or description, or the user can search for a table by entering a portion of the table description in the search box. Upon choosing a lookup table, the user will then have the option to edit, delete, or add records to the table.





Figure 53 CARIM Utilities

Please note that the LuTables (Look Up tables) perform a different function. Adding or editing the records in these tables will affect the data displayed within the system. Edits made to the table names will only affect the Utilities Module, and will not remove or change the table names in the database.


System Description CARIM

# **Appendix – Functional Process Flowcharts**


# CARIM High Level System Flowchart




# Receiving Module



# **Shipping Module**



# Inventory Module



# Entry from Main Menu

# **Utility Table Information Module**

