

Potency Assays for Existing Vaccines

John P. Hennessey, Jr.


Bioprocess Research & Development

Merck & Co., Inc

Assaying Potency of Novel Vaccines

Oct. 11 – 12, 2005

Potency: what are we after?



A Limited Survey of Licensed Vaccines

- US 52 vaccines licensed (CBER web site)
at least 4 more currently under review
28 product types
3 new types under review
- EU 13 vaccines centrally licensed (EMA web site)
at least 3 more currently under review
28 product types licensed in UK (MHRA web site)
3 new types under review

A Limited Survey of Licensed Vaccines

- Licensed (and under review) vaccines cover 24 different pathogens (59 if all serotypes counted)
 - 14 viral (24 types)
 - Most attenuated; some inactivated/recombinant
 - 10 bacterial (35 types)
 - Most purified components; some attenuated

A Limited Survey of Licensed Vaccines

- Licensed (and under review) vaccines cover 24 different pathogens (59 if all serotypes counted)

- 14 viral (24 types)

- | | | |
|--------------------------|---|--|
| • measles | influenza [#] (2) | hepatitis A [#] |
| • mumps | rabies [#] | hepatitis B [#] |
| • rubella | varicella | Japanese encephalitis virus [#] |
| • polio [#] (3) | small pox [#] | yellow fever virus [#] |
| • rotavirus (5)* | human papilloma virus [#] (4)* | |

- 10 bacterial (35 types)

* Under review


inactivated or recombinant form available

A Limited Survey of Licensed Vaccines


- Licensed (and under review) vaccines cover 24 different pathogens (59 if all serotypes counted)
 - 14 viral (24 types)
 - 10 bacterial (35 types)
 - *C. diphtheriae* *C. tetani* *B. pertussis*
 - *N. meningitidis* (4) *M. bovis** *S. typhii**
 - *S. pneumoniae* (23) *H. influenzae* B *V. cholerae*
 - *B. anthracis*

* attenuated forms available

Dissection of Vaccine Targets




Type of Potency Assay Used Reflects Product Type




Replication Assays

- Reflects the *in vitro* replication of the viable components (virus or bacteria) in a cell culture system
 - **Virus replication** typically reflects cytopathic effect (CPE) on (*i.e.* death of) cell substrate
 - Other approaches possible (to follow)
 - **Bacterial replication** typically reflects formation of visible colonies on a solid media substrate


In vitro Potency of Live Virus Vaccines



Virus Replication: Quantitative Plaque Formation Assay



Varicella Plaque Assay: Infected plate




Benefits:

- Simple
- Easily transferable procedure


Drawbacks:

- Limited specificity
- Manually intensive
- Poor precision/reproducibility
- Limited throughput

Virus Replication: TCID₅₀ Assay




- Benefits
 - Semi-automated
- Drawbacks
 - Limited specificity
 - Wide virus dilution ranges needed
 - Poor precision/reproducibility
 - Moderate throughput



Readout: TCID₅₀/mL

Root Assay Variability: ~60% RSD

Virus Replication: Quantitative PCR-based Potency Assay (QPA)



- Benefits:
 - Direct detection of viral (genome) replication
 - High throughput
 - Excellent specificity & sensitivity
 - Improved precision
- Drawbacks:
 - Costly system & reagents
 - New technology for QC lab

Bacterial Replication: Quantitative Colony Formation Assay

- Typically applying serial dilutions to solid-phase media, incubating, counting colonies
 - Bacterial equivalent of virus plaque assay
 - Specificity enhanced by use of selective media and evaluation of colony morphology
- Benefits and drawbacks similar to virus plaque assays

Immunogenicity Assay

- Reflects the *in vivo* biological response to the product in a relevant animal model
 - “Relevance” is relative
- Biomarker model
 - Assess biomarker (e.g. antibody levels) in pre- and post-vaccination sera
 - e.g. pertussis, Hib conjugate (no longer used)
 - Assess biomarker relative to reference standard
 - e.g. polio
- Challenge model
 - Neutralizing activity of serum from vaccinated animals in challenged animals
 - e.g. Diphtheria, tetanus

Immunogenicity Assay

- Benefits
 - Shows *in vivo* response
- Drawbacks
 - Relevance to predicting human response?
 - Requires continual supply of animals
 - Labor intensive
 - High variability


Antigenicity Assay

- Can be applicable to almost all vaccine components
 - Typically not used for “live” products
- Several formats to consider
 - Typical ELISA; extrapolation from std. curve
 - e.g. hepatitis A
 - IVRP; parallel line analysis, relative potency
 - e.g. hepatitis B, HPV
 - Competitive; four-parameter fit; relative potency
 - e.g. HPV (secondary assay)
 - Others (e.g. rate neph., rocket IE)
 - e.g. pneumococcal polysaccharides

Antigenicity Assay

- Benefits
 - Reflection of epitope content/integrity
 - Automatable
 - Relatively precise, moderate throughput
- Drawbacks
 - Presumed correlate of immunogenic potential
 - Sometimes bridged to immunogenicity assay (*e.g.* HPV)
 - Relevance to predicting human response?
 - Identifying/supplying reagents can be difficult


Structural Model of an HPV VLP



(Crystal structure coordinates courtesy of Prof. S. C. Harrison, Harvard University)


High-Resolution AFM Image of an HPV VLPs

HPV 16




Morley & Y. Wang

HPV 18




(20-Å Resolution *False-color* Image, Cryo-TEM)
Watson (AP-MSD), Duncan (MRL) & Potter (Scripps)

Correlation of Mouse Immunogenicity with Relative Antigenicity (IVRP)



Correlation of Potency Measures for All Constructs



Neutralization Assay

- Multiple approaches to “neutralization”
- Use specific antisera to block response
 - Replication assay
 - “Neutralize” replication
 - typically used for “identity” rather than potency
 - Immunogenicity assay
 - Immune response “neutralizes” pathogenic response
 - e.g. animal challenge assay (diphtheria, tetanus)
 - Antigenicity assay
 - “Neutralize” epitope binding
 - typically used for “identity” rather than potency


Mass Assay

- Assay that reflects “mass” of component in product (assumes potency-to-mass ratio is fixed)
 - Colorimetric assay
 - polysaccharide content (*S. typhii* Ps)
 - Secondary evaluation of MW and O-Ac content
 - protein assay (hepatitis B; secondary assay)
 - HPLC assay
 - polysaccharide content (Hib conjugate)
 - Secondary evaluation of MW and % conjugated

Mass Assay

- Benefits
 - Physico-chemical measurement
 - Limited variability
 - Automatable
 - Relatively precise, moderate throughput
- Drawbacks
 - High burden of proof that mass \sim potency
 - Relevance to predicting human response?

Choice of Potency Assay Depends on Product Type



New vaccine types (e.g. viral vectors, plasmids) will require consideration of new potency assay types.